

Titration and measuring

Method Alkalinity_Op
Method saving date 2018-06-06 14:41:49 UTC+3
Method version 2
Method group Main group
Method status original
Method saved by (full name)
Method saved by (short name) Owner

MEAS pH

MEAS init pH

General/Hardware

Device

Device name 855_1
Device type 855 Robotic Titrosampler

Sensor

Measuring input 1
Sensor pH electrode
Temperature measurement automatic

Stirrer

Stirrer 1
Stirring rate 5
Switch off automatically on

Measuring parameters

Measurement

Measurement with drift control on
Signal drift 1 mV/min
Min. waiting time 30 s
Max. waiting time 154 s
Measuring interval 2.0 s
Stop measured value pH off
Measurement without drift control off

Temperature

Temperature 20 °C

Evaluations

Fixed endpoint evaluation off
Minimum evaluation off
Maximum evaluation off
Break point evaluation off

Additional measured values

Additional calculated measured values off
Additional external measured values off

DET pH

DET pH

General/Hardware

Device

Device name 855_1

Device type 855 Robotic Titrosampler

Dosing device

Dosing device 2

Solution HCl

Sensor

Measuring input 1

Sensor pH electrode

Temperature measurement automatic

Stirrer

Stirrer 1

Stirring rate 5

Switch off automatically on

Start conditions

Initial measured value

Signal drift off mV/min

Min. waiting time 0 s

Max. waiting time 1 s

Start volume

Start volume 1.5 mL

Dosing rate maximum mL/min

Start measured value

Start measured value pH off

Dosing rate 5 mL/min

Start slope

Start slope off pH/mL

Dosing rate 5 mL/min

Pause

Pause 0 s

Titration parameters

Titration rate

Titration rate fast

Temperature

Temperature 20 °C

Stop conditions

Stop volume 10 mL

Stop measured value pH 2.8

Stop EP 1

Volume after EP 1 mL

Stop time off s

Filling rate maximum mL/min

Potentiometric evaluation

Evaluation without window off

Evaluation with measured value window (pH) on

Lower limit pH	Upper limit pH	EP criterion	EP recognition
3	5	7	greatest

Evaluation with volume window (mL) off

Additional evaluations

Fixed endpoint evaluation off

pK/HNP evaluation off

Minimum evaluation off

Maximum evaluation off

Break point evaluation off

Gran evaluation on

Procedure Standard

Initial volume 0.2 mL

Lower limit pH -20.0

Upper limit pH 20.0

Additional measured values

Additional calculated measured values off

Additional external measured values off