

Method parameters

Method Alkalinity_Op
Method saving date 2018-06-06 14:41:49 UTC+3
Method version 2
Method group Main group
Method status original
Method saved by (full name)
Method saved by (short name) Owner

START

Main track

General

Workplace view
Current view on
Track view for live window
Live display 1 Main track
Live display 2 Main track
Statistics off
Number of single determinations not defined
Conditioning
Automatic conditioning on

Application note

Method variables

Name	Type	Assignment	Fixed value	Comment	Monitoring
Sample size	Number	Sample size		Sample size	off
Sample size unit	Text	Sample size unit		Sample size unit	off
Sample position	Number	Sample position		Sample position number	off
ID1	Text	ID1		Sample identification 1	off
ID2	Text	ID2		Sample identification 2	off
ID3	Text	ID3		Sample identification 3	off

Name **Sample size**
Type Number
Assignment on. Sample size
Fixed value off.
Check at start on
Comment Sample size
Variable monitoring off
Lower limit
Upper limit
Message
Display message on
Record message on

Message by e-mail off
 Use e-mail template off
 E-mail template
 Mail to
 Subject Message from tiamo - Method '%1' - Command '%2'
 User
 Mail from
 SMTP Server
 POP3 Server
 Acoustic signal off
 Action off
 Cancel determination on
 Cancel determination and series off

Name **Sample position**
 Type Number
 Assignment on. Sample position
 Fixed value off.
 Check at start on
 Comment Sample position number
 Variable monitoring off
 Lower limit
 Upper limit
 Message
 Display message on
 Record message on
 Message by e-mail off
 Use e-mail template off
 E-mail template
 Mail to
 Subject Message from tiamo - Method '%1' - Command '%2'
 User
 Mail from
 SMTP Server
 POP3 Server
 Acoustic signal off
 Action off
 Cancel determination on
 Cancel determination and series off

Name **Sample size unit**
 Type Text
 Assignment on. Sample size unit
 Fixed value off.
 Check at start on

Comment Sample size unit

Name **ID1**
 Type Text
 Assignment on. ID1
 Fixed value off.
 Check at start on
 Comment Sample identification 1

Name **ID2**
 Type Text
 Assignment on. ID2
 Fixed value off.
 Check at start on
 Comment Sample identification 2

Name **ID3**
 Type Text
 Assignment on. ID3
 Fixed value off.
 Check at start on
 Comment Sample identification 3

SEQUENCE SEQUENCE 5

MOVE

MOVE 6

Device
 Device name 855_1
 Device type 855 Robotic Titrosampler
 Target
 Tower 1
 Move Sample position
 Beaker test
 Display message on
 Cancel determination off
 Cancel determination and series off
 Parameters
 Shift rate 20 %/s
 Shift direction auto
 Swing rate 55 %/s

LIFT

LIFT 7

Device
 Device name 855_1
 Device type 855 Robotic Titrosampler
 Target
 Tower 1

Lift position Work position mm
Parameters
Lift rate 25 mm/s

MEAS pH**MEAS init pH****General/Hardware**

Device

Device name 855_1
Device type 855 Robotic Titrosampler

Sensor

Measuring input 1
Sensor pH electrode
Temperature measurement automatic

Stirrer

Stirrer 1
Stirring rate 5
Switch off automatically on

Measuring parameters

Measurement

Measurement with drift control on
Signal drift 1 mV/min
Min. waiting time 30 s
Max. waiting time 154 s
Measuring interval 2.0 s
Stop measured value pH off
Measurement without drift control off

Temperature

Temperature 20 °C

Evaluations

Fixed endpoint evaluation off
Minimum evaluation off
Maximum evaluation off
Break point evaluation off

Additional measured values

Additional calculated measured values off
Additional external measured values off

DET pH**DET pH****General/Hardware**

Device

Device name 855_1
Device type 855 Robotic Titrosampler

Dosing device

Dosing device 2
Solution HCl

Sensor

Measuring input 1
Sensor pH electrode
Temperature measurement automatic
Stirrer
Stirrer 1
Stirring rate 5
Switch off automatically on

Start conditions

Initial measured value
Signal drift off mV/min
Min. waiting time 0 s
Max. waiting time 1 s
Start volume
Start volume 1.5 mL
Dosing rate maximum mL/min
Start measured value
Start measured value pH off
Dosing rate 5 mL/min
Start slope
Start slope off pH/mL
Dosing rate 5 mL/min
Pause
Pause 0 s

Titration parameters

Titration rate
Titration rate fast
Temperature
Temperature 20 °C

Stop conditions

Stop volume 10 mL
Stop measured value pH 2.8
Stop EP 1
Volume after EP 1 mL
Stop time off s
Filling rate maximum mL/min

Potentiometric evaluation

Evaluation without window off
Evaluation with measured value window (pH) on

Lower limit pH	Upper limit pH	EP criterion	EP recognition
3	5	7	greatest

Evaluation with volume window (mL) off

Additional evaluations

Fixed endpoint evaluation off

pK/HNP evaluation off
Minimum evaluation off
Maximum evaluation off
Break point evaluation off
Gran evaluation on
 Procedure Standard
 Initial volume 0.2 mL
 Lower limit pH -20.0
 Upper limit pH 20.0

Additional measured values

Additional calculated measured values off
Additional external measured values off

SEQUENCE SEQUENCE 8

MOVE

MOVE 11

Device

Device name 855_1
Device type 855 Robotic Titrosampler

Target

Tower 1
Move Special beaker
Number 1

Beaker test

Display message on
Cancel determination off
Cancel determination and series off

Parameters

Shift rate 20 °/s
Shift direction auto
Swing rate 55 °/s

LIFT

LIFT 12

Device

Device name 855_1
Device type 855 Robotic Titrosampler

Target

Tower 1
Lift position Rinse position mm

Parameters

Lift rate 25 mm/s

PUMP

PUMP 15

Device

Device name 855_1
Device type 855 Robotic Titrosampler

Pumps

Tower 1
Pump(s) 1
Action
Switch on off
Switch off off
Duration on
Time 3 s

EXIT **Exit track**

CALC **CALC**

Result name	Formula	Unit	Decimal places	Assignment	Statistics
Alkalinity	$=(\text{'DET pH.EP\{1\}.VOL'}/1000) * \text{'DET pH.CONC'*1000000 / ('MV.Sample size'}/1000)$	μmol/kg	3	RS01	off
dKH	$= \text{'RS.Alkalinity'} * 0.00205 / 2 * 2.8$	dKH	2	RS05	off
pH at EP	$= \text{'DET pH.EP\{1\}.MEA'}$		2	RS02	off
pH init	$= \text{'MEAS init pH.EME'}$		2	RS03	off
pH final	$= \text{'DET pH.EME'}$		2	RS04	off
Alkalinity_2	$= \text{'RS.Alkalinity'} * 0.00205 / 2 * 2.8*1.025*0.357$	Meq/kg	3	RS06	off

Result name **Alkalinity**
Formula $=(\text{'DET pH.EP\{1\}.VOL'}/1000) * \text{'DET pH.CONC'*1000000 / ('MV.Sample size'}/1000)$
Unit μmol/kg
Decimal places 3
Assignment RS01
Statistics off
Description
Result monitoring off
Save result as common variable off
Common variable
Save result as titer off
Solution name HCl

Result name **dKH**
Formula $= \text{'RS.Alkalinity'} * 0.00205 / 2 * 2.8$
Unit dKH
Decimal places 2
Assignment RS05
Statistics off
Description
Result monitoring off

Save result as common variable off
Common variable
Save result as titer off
Solution name HCl

Result name **pH at EP**
Formula = 'DET pH.EP{1}.MEA'
Unit
Decimal places 2
Assignment RS02
Statistics off
Description
Result monitoring off
Save result as common variable off
Common variable
Save result as titer off
Solution name HCl

Result name **pH init**
Formula = 'MEAS init pH.EME'
Unit
Decimal places 2
Assignment RS03
Statistics off
Description
Result monitoring off
Save result as common variable off
Common variable
Save result as titer off
Solution name HCl

Result name **pH final**
Formula = 'DET pH.EME'
Unit
Decimal places 2
Assignment RS04
Statistics off
Description
Result monitoring off
Save result as common variable off
Common variable
Save result as titer off
Solution name HCl

Result name **Alkalinity_2**
Formula = 'RS.Alkalinity' * 0.00205 / 2 * 2.8*1.025*0.357

Unit Meq\kg
Decimal places 3
Assignment RS06
Statistics off
Description
Result monitoring off
Save result as common variable off
 Common variable
Save result as titer off
 Solution name HCl

REPORT

REPORT

Report template
 Report template example report
Report output
 Printer off
 PDF file off

DATABASE

DATABASE

Database

tiamo

EXPORT

EXPORT 18

Export template
 Export template O. patagonica physiology

SERIES START

Series start track

RACK

initialize rack

Device
 Device name 855_1
 Device type 855 Robotic Titrosampler
Rack test off

MOVE

to rinse beaker

Device
 Device name 855_1
 Device type 855 Robotic Titrosampler
Target
 Tower 1
 Move Special beaker
 Number 1
Beaker test
 Display message off
 Cancel determination on

Cancel determination and series off

Parameters

Shift rate 20 %/s

Shift direction auto

Swing rate 55 %/s

LIFT

LIFT 2

Device

Device name 855_1

Device type 855 Robotic Titrosampler

Target

Tower 1

Lift position Rinse position mm

Parameters

Lift rate 25 mm/s

PUMP

aspirate and rinse

Device

Device name 855_1

Device type 855 Robotic Titrosampler

Pumps

Tower 1

Pump(s) 1

Action

Switch on off

Switch off off

Duration on

Time 7 s

**SERIES
END**

Series end track

MOVE

MOVE 16

Device

Device name 855_1

Device type 855 Robotic Titrosampler

Target

Tower 1

Move Special beaker

Number 2

Beaker test

Display message on

Cancel determination off

Cancel determination and series off

Parameters

Shift rate 20 %/s

Shift direction auto

Swing rate 55 %/s

LIFT**LIFT 1**

Device

Device name 855_1

Device type 855 Robotic Titrosampler

Target

Tower 1

Lift position Rinse position mm

Parameters

Lift rate 25 mm/s

STIR**STIR 17**

Device

Device name 855_1

Device type 855 Robotic Titrosampler

Stirrer

Stirrer 1

Stirrer type unknown

Stirring rate 3

Action

Switch on off

Switch off off

Duration on

Time 30 s